

AFIAS Myoglobin

INTENDED USE

AFIAS Myoglobin is a fluorescence immunoassay (FIA) for the quantitative determination of Myoglobin in human whole blood/serum/plasma. It is useful as an aid in management and monitoring of acute myocardial infarction (AMI).
For *in vitro* diagnostic use only.

INTRODUCTION

Myoglobin is an iron- and oxygen-binding protein found in both skeletal and myocardial muscles. It acts as a transport protein and is involved in diffusion of oxygen in the muscle tissue. Myoglobin is a single-chain globular protein of 154 amino acids. It is composed of a central iron-containing ‘Heme’ which is enclosed in a compact bundle-like or prism-like arrangement formed by the eight right-handed α -helices^{1,2}. Being a cytoplasmic protein having low molecular weight (of 17,699 daltons), myoglobin is released into the serum more rapidly as compared to other cardiac markers upon damage to the myocardial cells. Serum concentration of myoglobin increases above the normal range as early as 1 hour after acute myocardial infarction (AMI), attains peak level in approximately 4 to 8 hours after the onset and normalize rapidly afterwards. Thus, myoglobin is better suited as a cardiac marker for early diagnosis of AMI. However, the elevated myoglobin is not specific to AMI owing to its large quantities in skeletal muscles as well. Despite its low clinical specificity and weak predictive value towards AMI, myoglobin is still a promising cardiac marker when other markers such as Creatin Kinase Isoenzyme-MB (CK-MB) and Cardiac Troponin-I (cTn-I) as well as other indicators like clinical signs and ECG are taken into account for diagnosis/confirmation of AMI³⁻⁸

PRINCIPLE

The test uses a sandwich immunodetection method.
The detector antibodies in buffer bind to antigens in the sample, forming antigen-antibody complexes, and migrate onto nitrocellulose matrix to be captured by the other immobilized antibodies on a test strip.
More antigens in the sample will form more antigen-antibody complexes which lead to stronger fluorescence signal by detector antibodies, which is processed by the instrument for AFIAS tests to show Myoglobin concentration in the sample.

COMPONENTS

AFIAS Myoglobin consists of ‘cartridges’.

- Each sealed aluminum pouch contains two cartridges.
- Each cartridge packaged in an aluminum pouch has three components including a cartridge part, a detector part and a diluent part.
- The cartridge part contains the membrane called a test strip which has anti-myoglobin at the test line, and streptavidin at the control line.
- The detector part contains anti-myoglobin-fluorescence

conjugate, biotin-BSA-fluorescence conjugate and sodium azide as a preservative in phosphate buffered saline (PBS).

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- Follow instructions and procedures described in this ‘Instructions for use’.
- Use only fresh samples and avoid direct sunlight.
- Lot numbers of all the test components (cartridge and ID chip) must match each other.
- Do not interchange the test components between different lots or use the test components after the expiration date, either of which might yield incorrect test result(s).
- Do not reuse cartridges. A cartridge should be used for testing one sample only.
- The cartridge should remain sealed in its original pouch until just before use. Do not use a cartridge, if the pouch is damaged or has already been opened.
- Frozen sample should be thawed only once. For shipping, samples must be packed in accordance with local regulations. Sample with severe hemolysis and/or hyperlipidemia must not be used.
- If test components and/or sample are stored in refrigerator, then allow cartridge and sample to be at room temperature for approximately 30 minutes before use.
- The instrument for AFIAS tests may generate slight vibration during use.
- The cartridge contains sodium azide (NaN₃), and it may cause certain health issues like convulsions, low blood pressure, low heart rate, loss of consciousness, lung injury and respiratory failure. Avoid contact with skin, eyes, and clothing. In case of contact, rinse immediately with running water.
- No Biotin interference was observed in **AFIAS Myoglobin** when biotin concentration in the sample was below 200 ng/mL. If a patient has been taking biotin at dosage of more than 0.03 mg a day, it is recommended to test again 24 hours after discontinuation of biotin intake.
- AFIAS Myoglobin** will provide accurate and reliable results subject to the below conditions.
 - AFIAS Myoglobin** should be used only in conjunction with the instrument for AFIAS tests.
 - Have to use recommended anticoagulant.

Recommended anticoagulant
K ₃ EDTA, Lithium heparin, Sodium citrate

LIMITATIONS OF THE TEST SYSTEM

- The test may yield false positive result(s) due to the cross-reactions and/or non-specific adhesion of certain sample components to the capture/detector antibodies.
- The test may yield false negative result(s) due to the non-responsiveness of the antigens to the antibodies which is the most common if the epitope is masked by some unknown components, so therefore not being able to be detected or captured by the antibodies. The instability or degradation of the antigens with time and/or temperature may also cause false negative result as it makes the antigens unrecognizable by the antibodies.
- Other factors may interfere with the test and cause erroneous results, such as technical/procedural errors, degradation of the test components/reagents or presence of interfering substances in the test samples.

- Any clinical diagnosis based on the test result must be supported by a comprehensive judgment of the concerned physician in conjunction with clinical symptoms and other relevant test results.

STORAGE AND STABILITY

Storage condition			
Component	Storage Temperature	Shelf life	Note
Cartridge	2 - 8 °C	20 months	Unopened
		1 month	Resealed

- Return an unused cartridge to the spare cartridge zipper bag containing the desiccant pack. Reseal along entire edge of zip-seal.

MATERIALS SUPPLIED

REF SMFP-34 Components of AFIAS Myoglobin	
Cartridge box:	
- Cartridge	24
- Pipette tip (zipper bag)	24
- Spare cartridge zipper bag	1
- ID chip	1
- Instructions for use	1

MATERIALS REQUIRED BUT SUPPLIED ON DEMAND

Following items can be purchased separately with **AFIAS Myoglobin**.
Please contact our sales division for more information.

Instrument for AFIAS tests	
- AFIAS-1	REF FPRR019
- AFIAS-3	REF FPRR040
- AFIAS-6	REF FPRR020
- AFIAS-10	REF FPRR038
Boditech Cardiac Control	REF CFPO-98
Boditech Cardiac Calibrator	REF CFPO-110

SAMPLE COLLECTION AND PROCESSING

- The sample type for **AFIAS Myoglobin** is human whole blood/serum/plasma.
- It is recommended to test the sample within 24 hours after collection when collected sample is stored at room temperature.
 - The samples (serum, plasma) should be separated from the clot by centrifugation within 3 hours after the collection of whole blood.
 - The samples (whole blood, serum, plasma) may be stored for a week at 2-8 °C prior to being tested. If testing will be delayed more than a week, samples (serum, plasma) should be frozen at -20 °C.
 - The samples (serum, plasma) stored frozen at -20 °C for 3 months showed no performance difference.
 - However, the whole blood sample should not be kept in a freezer in any case.
 - As a repeated freeze-thaw cycle may affect the test result, do not refreeze previously frozen samples.

TEST SETUP

- Check the components of the **AFIAS Myoglobin** as described below. : Cartridges, pipette tips, an ID chip, a spare cartridge zipper bag and an instructions for use.
- Ensure that the lot number of the cartridge matches that of the ID chip.
- If the sealed cartridge has been stored in a refrigerator, place them on a clean and flat surface at room temperature for at least 30 minutes before testing.
- Turn on the instrument for AFIAS tests.
- Empty the tip box.
- Insert the ID chip into the ‘ID chip port’.

※ **Please refer to the instrument for AFIAS tests operation manual for complete information and operating instructions.**

TEST PROCEDURE

► AFIAS-1, AFIAS-3, AFIAS-6

General mode

- Insert a cartridge into the cartridge holder.
- Insert a tip into the tip hole of the cartridge.
- Select the ‘General mode’ in the instrument for AFIAS tests.
- Take 100 μ L of sample (whole blood/serum/plasma/control) using a pipette and dispense it into the sample well of the cartridge.
- Tap the ‘Start’ button on the screen.
- The test result will be displayed on the screen after 12 minutes.

► AFIAS-10

Normal mode

- Insert a cartridge into the cartridge holder.
- Insert a tip into the tip hole of the cartridge.
- Tap the ‘Load’ button of the bay that holds the cartridge with the tip to read the barcode of the cartridge and please confirm the item name written on the cartridge.
- Insert the sample tube into the tube rack.
- Insert the tube rack into the loading part of the sampling station.
- Tap the ‘Start’ button on the screen.
- The test result will be displayed on the screen after 12 minutes.

Emergency mode – General tip

- The test procedure is same with the ‘Normal mode 1) – 3)’.
- Convert the ‘Emergency mode’ in AFIAS-10.
- Select the tip type (general tip) on the screen.
- Select the sample type (whole blood/serum/plasma) on the screen.
- Take 100 μ L of the sample using a pipette and dispense it into the sample well of the cartridge.
- Tap the ‘Start’ button on the screen.
- The test result will be displayed on the screen after 12 minutes.

INTERPRETATION OF TEST RESULT

- The instrument for AFIAS tests calculates the test result automatically and displays myoglobin concentration of the test sample in terms of ng/mL.
- Reference value: 70 ng/mL
- Working range: 2-500 ng/mL

QUALITY CONTROL

- Quality control tests are a part of the good testing practice to confirm the expected results and validity of the assay and should be performed at regular intervals.
- Quality control tests should also be performed whenever there is any question concerning the validity of the test results.
- Control materials are provided on demand with **AFIAS Myoglobin**. For more information regarding obtaining the control materials, contact Boditech Med Inc.'s Sales Division for assistance.
(Please refer to the instructions for use of control material.)

PERFORMANCE CHARACTERISTICS

- **Analytical sensitivity**
Limit of Blank (LoB) 1.30 ng/mL
Limit of Detection (LoD) 1.95 ng/mL
Limit of Quantitation (LoQ) 2.00 ng/mL
- **Analytical specificity**
 - Cross-reactivity
Biomolecules such as below the ones in the table were added to the test sample(s) at concentrations much higher than their normal physiological levels in the blood. **AFIAS Myoglobin** test results did not show any significant cross-reactivity with these biomolecules.
- Interference
Interferents listed in the following table were added to the test sample at the concentration mentioned below. **AFIAS Myoglobin** test results did not show any significant interference with these materials.

Cross-reactants	Concentration
Troponin complex	1,000 ng/mL
CK-MB	1,000 ng/mL
D-Dimer	1,000 ng/mL
NT-proBNP	1,000 ng/mL

Interferents	Concentration
D-glucose	55.5 mmol/L
L-Ascorbic acid	175 µmol/L
Bilirubin (unconjugated)	684 µmol/L
Hemoglobin	10 g/L
Cholesterol	10.3 mmol/L
Triglyceride	16.94 mmol/L
Heparin	330 U/dL
EDTA	3.4 µmol/L
sodium citrate	2 mg/mL

Precision

- Single-site study
Repeatability (within-run precision)
within-laboratory precision (Total precision)
Lot to lot precision
3 Lots of **AFIAS Myoglobin** were tested for 20 days. Each standard material was tested 2 times per day. For each test, each material was duplicated.

Myoglobin [ng/mL]	Single-site study					
	Repeatability		within-laboratory precision		Lot to lot precision	
	AVG [ng/mL]	CV (%)	AVG [ng/mL]	CV (%)	AVG [ng/mL]	CV (%)
55	55.03	6.18	55.14	5.91	55.05	5.86
100	98.91	6.21	100.25	6.14	99.68	5.88
300	297.45	5.73	298.69	5.40	299.20	5.57

- Multi-site study
Reproducibility
1 Lot of **AFIAS Myoglobin** was tested for 5 days in 3 different sites (1 person per 1 site, 1 instrument per 1 site). Each standard material was tested 1 time per and 5 replicates per day.

Myoglobin [ng/mL]	Multi-site study	
	AVG [ng/mL]	CV (%)
55	54.35	6.2
100	100.36	6.0
300	302.48	5.9


Accuracy

- The accuracy was confirmed by testing with 3 different lots of **AFIAS Myoglobin**. The tests were repeated 10 times at each concentration of the control standard.

Myoglobin [ng/mL]	LOT 1	LOT 2	LOT 3	AVG [ng/mL]	Recovery (%)
2.00	2.08	2.09	2.07	2.08	104
101.60	101.61	104.37	103.54	103.17	102
201.20	206.26	192.36	200.63	199.75	99
300.80	308.92	304.16	306.06	306.38	102
400.40	405.11	396.91	397.00	399.67	100
500.00	487.28	480.77	477.27	481.77	96

Comparability

Myoglobin concentration of 120 clinical samples were independently with **AFIAS Myoglobin (AFIAS-6)** and **comparator A** as per prescribed test procedures. Test results were compared, and their comparability was investigated with linear regression and correlation coefficient (R). The regression equation and correlation coefficient are as follow.



REFERENCES

1. Kent Lewandrowski, Ahchean Chen and James Januzzi, Cardiac markers for myocardial infarction, Am J Clin Pathol 2002; 118 (Suppl 1):S93-S99.
2. Cox, MM, Nelson, DL. Lehninger: Principles of Biochemistry, 3rd edition. W.H. Freeman and Company, New York, 2000, 206.
3. Ordway GA, Garry DJ. Myoglobin: An essential hemoprotein in striated muscle. J Exp Biol.. 2004; 207(Pt 20):3441-6.
4. Lewandrowski K, Chen A, Januzzi J. Cardiac markers for myocardial infarction. A brief review. Am J Clin Pathol. 2002;118:S93-9.
5. Vaidya HC. Myoglobin: an early biochemical marker for the diagnosis of acute myocardial infarction. J Clin Immunoassay. 1994; 17:35-39.
6. Gibler WB, Gibler CD, Weinshenker C, et al. Myoglobin as an early indicator of acute myocardial infarction. Ann Emerg Med. 1987; 16:851-856.
7. Adams JE, Abendschein DR, Jaffe AS. Biochemical markers of myocardial injury: is MB creatine kinase the choice for the 1990s? Circulation. 1993; 88:750-763.
8. Mair J, Morandell D, Genser N, et al. Equivalent early sensitivities of myoglobin, creatine kinase-MB mass, creatine kinase isoforms ratios, and cardiac troponins I and T for acute myocardial infarction. Clin Chem. 1995; 41:1266-1272.
9. Mercer DW. Role of cardiac markers in evaluation of suspected myocardial infarction. Postgrad Med. 1997;102:113-122

Note: Please refer to the table below to identify various symbols.

	Sufficient for <n> tests
	Read instruction for use
	Use by Date
	Batch code
	Catalog number
	Caution
	Manufacturer
	Authorized representative of the European Community
	In vitro diagnostic medical device
	Temperature limit
	Do not reuse
	This product fulfills the requirements of the Directive 98/79/EC on in vitro diagnostic medical devices

For technical assistance, please contact:

Boditech Med Inc.'s Technical Services

Tel: +(82) -33-243-1400

E-mail: TS@boditech.co.kr

Boditech Med Inc.

43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gang-won-do, 24398, Republic of Korea

Tel: +(82) -33-243-1400

Fax: +(82) -33-243-9373

www.boditech.co.kr

Obelis s.a

Bd. Général Wahis 53, 1030 Brussels, Belgium

Tel: +(32) -2-732-59-54

Fax: +(32) -2-732-60-03

E-Mail: mail@obelis.net