IgG id C€

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REF

3300505

IgG

1 x 15 tests

Sólo para uso diagnóstico *in vitro*

IgG

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Mancini et al., Immunochemistry, 2: 235 (1965)

IgG id

CONTENTS			
REF	3300505	IgG	1 x 15 tests
For <i>in vitro</i> diagnostic use only			

IgG

Quantitative determination of IgG by radial immunodiffusion

PRINCIPLE

The examined protein, diffusing in agarose gel containing a specific antibody will form an immuno-complex, visible as a ring around the well. The ring diameter is in proportion to the concentration of the analyzed protein. The proportion corresponds to the diffusion time. In fact, at the end (72 h), the square of the diameter will be in linear proportion to the concentration (*Procedure 1*), while after a shorter period of time the square of the diameter will be in a logarithmic proportion to the concentration (*Procedure 2*). In both cases, a calibration curve should be constructed, using at least three calibration points. However, a reference table is provided showing the relation between any concentration and the diameter at the end of the analysis time.

REAGENT COMPOSITION

Plate with agarose gel containing the antiserum for the specific protein to be analyzed.

STORAGE AND STABILITY

 Store at 2-8°C.
The plates are ready to be used.
The Reagents are stable until the expiry date stated on the label.

SAMPLES

Serum or plasma. Stable six days at 2-8 °C.

ASSAY PROCEDURE

Allow the plate to come to room temperature, open it and, if moisture is present on the wells, wait until it has evaporated. Withdraw 5 µL of samples and/or controls and wait until it has been completely absorbed before handling the plate. Cover it tightly and place it in a moist chamber during the incubation time:

- 72 hours for the *Procedure 1* or 3.
- 18 hours for the *Procedure 2*.

READING TEST RESULTS

Measure the precipitation ring to the nearest 0.1 mm, after the required period according to used assay procedure, 72 h for *Procedure 1* or 3, 18 h for *Procedure 2*.

Procedure 1

Reproduce on a plotting paper, the square of the precipitation rings and the concentration of the controls. A line with an intercept should be in the range 10-12 mm². Sample values are determined by interpolation.

Procedure 2

Reproduce on a plotting paper, the square of the precipitation rings against their concentration logarithms. The graph will be a straight line only for a low values. For high values may curve a little. Sample values are obtained by interpolation.

Procedure 3

Read on the enclosed Reference Table, the concentration value corresponding to the precipitation ring diameter. The ring value obtained for the control should have a confidence limit of 0.2 mm (from the reference table).

NOTES

The diffusion time and the reading time depend on the concentration and the specific diffusing protein. After 72 hours, the diffusion of the protein, at any concentration, is completed. For lower concentrations, it is possible to read at lower times (ie 36 h). However, in such cases it is advisable to read again after 3-5 hours. If the diameter is still the same, it is possible to set the concentration, on the contrary, if the diameter is different, the ring should be remeasured after a further 3-5 hours.

QUALITY CONTROL

It is suggested to perform an internal quality control. Not included in the kit.

ANALYTICAL PERFORMANCE

- **Measure's limit.** 300-3500 mg/dL.
- **Precision**

mg/dL	Intra-assay		Inter-assay	
Mean	1299.80	1869.18	1302.32	1863.52
SD	20.78	26.68	21.90	25.23
CV%	1.60	1.43	1.68	1.35
N	10	10	20	20

Replicates: 10 for each level.

Replicates: 20 for each level.

- **Correlation.** This assay (x) was compared with a similar commercial method (y). The results were:

$$N = 70 \quad r = 0.99942 \quad y = 0.9835x + 23.161$$

REFERENCE VALUES

IgG 800-1800 mg/dL

REFERENCES

Fahey et al., *J. Immunol.*, 94: 84 (1965)
Mancini et al., *Immunochemistry*, 2: 235 (1965)

RI33005-2/1104